



Monitoring and Measuring Plate Flatness During Lapping



Lapping and
Polishing

1.0: Purpose

Lapping and polishing processes are typically performed on a hard, metal plate used in conjunction with abrasive suspensions such as diamond, silicon carbide (SiC), aluminum oxide (Al₂O₃), or boron carbide (B₄C). The metal lapping plate selected depends upon the desired material removal rate, the surface finish desired, the hardness of the workpiece being lapped, and the flatness requirement. During the lapping process material is removed from both the workpiece and the lapping plate, which eventually can lead to a lapping plate that is out of flatness. Understanding the dynamics of plate wear and the tendencies a plate has with respect to flatness can be helpful in process monitoring and general instrument maintenance. This paper describes a measurement technique used for monitoring the plate flatness at various stages to understand the mechanisms behind plate wear.

2.0: Experiments and Procedures

2.1: Measurement Technique

Establishing a measurement technique to evaluate the surface of a cast iron lapping plate, Model LP 920M, was first developed. The lapping plate was placed onto a three point adjustable table on the top of a granite leveling plate. A digital depth gauge (DDG) with 0.001" resolution was used in conjunction with an analog dial gauge (ADG) with 0.0001" resolution. The gauge was used to measure the plate surface at the three contact points. Based upon these three points the plate could be adjusted so that it was parallel. The gauge was then used to measure various points on the lapping plate surface and the data points were plotted to show total variation in plate flatness. This setup is illustrated schematically in Figure 1 below.




Figure 1: Schematic illustration of the technique used for measuring plate flatness. The plate is first made parallel on a granite leveling block using a dial gauge and three adjustment legs. After the plate is level, measurements are taken at various points on the surface of the lapping plate with the dial gauge. Each data point is then plotted to provide a surface plot illustrating the flatness of the plate after any given process.

2.3.1: Initial Plate Measurements




Figure 3: Plot of a new cast iron lapping plate prior to use. The plot shows three plotted series which represent the data series A, B, and C shown in Figure 1. From this plot you can see that the lapping plate is far from the tolerance specifications given for new lapping plates. The plate is dipped in the center, creating a concave plate shape.


From the plot above it can be easily seen that the lapping plate started off being concave in shape, exhibiting a severe flatness issue with the existing plate. These plates must be within the specified plate tolerances, otherwise poor specimen quality and unflat specimens will result.

2.3.1: 'Standard' Plate Conditioning Process

Plate conditioning was done using the 'standard' conditioning ring arrangement. The following parameters were used for conditioning:

- LAP SPEED: 2.5 (125 rpm)
- ARM SPEED: 10
- ABRASIVE: 23 μm BORON CARBIDE (B_4C)
- LAPPING PLATE: LP 920M FLAT CAST Fe
- CONDITIONING TIME: 30 MINUTES

CONDITIONING RING LOCATION







Figure 4: Plot of a new cast iron lapping plate following conditioning using the 'standard' conditioning ring method. As can be seen, the plate has flattened at the bottom, however no appreciable change in the plate shape can be observed. The plate still has a severe concave shape that will hinder proper specimen preparation.

2.3.2: Modified Plate Conditioning Process

Using a modified plate conditioning arrangement the lapping plate was conditioned for 30 minutes using the same parameters used in Section 2.3.1.

LAP SPEED: 2.5 (125 rpm)
 ARM SPEED: 10
 ABRASIVE: 23 μ m BORON CARBIDE (B₄C)
 LAPPING PLATE: LP 920M FLAT CAST Fe
 CONDITIONING TIME: 32 MINUTES

CONDITIONING RING LOCATION







Figure 5,6: Plot of cast iron lapping plate following conditioning with the modified arrangement. As can be seen in Figure 5 (above) the lapping plate has become close to the tolerance limits specified in the manufacturing drawings of the plate, indicating conditioning using the modified arrangement is much more effective in producing flat plates. In Figure 6 (below) the plate has been conditioned to be within the specifications of the lapping plate.



3.0: Lapping Experiments

Once a flat lapping plate had been produced, the effects of lapping a specimen using a Model 155 Lapping and Polishing Fixture were examined to determine the rate of material removal on the cast iron lapping plate during specimen processing. Understanding how long specimens take to make the plate unflat will help in determining the proper amount of time before plate conditioning is required.

A Model 155 Lapping and Polishing Fixture was set up with a series of glass/Si/glass sandwich samples for lapping. A total of four specimens were mounted for vertical edge polishing, with each specimen assembly approximately 10mm in thickness. Lapping was completed in 40 minutes using a 23 μ m boron carbide abrasive slurry. The load selected for lapping was 300 grams and the lap speed and arm speed were 2 and 10, respectively. Below is a schematic illustration of the lapping fixture position on the lapping plate (A), the plot of the surface (B), and the resulting plate surface after lapping was completed (C).




Figure 6: Illustration of the lapping arrangement when using the Model 155 Lapping and Polishing Fixture on the Model 920. In A) the Model 155 is being used as normal. In B) the resulting plate surface plot after lapping with the lapping fixture. In C) a schematic illustration of the resulting plate shape after lapping. The lapping fixture removes material faster where it is being used and the areas not exposed to the fixture are not removed as quickly. Although the plate is still relatively flat, conditioning is needed to get the plate back into its specified range of flatness.

3.1: Plate Re-conditioning

Following lapping the plate was re-conditioned using the modified conditioning ring arrangement. Tracking the behavior of the plate was done by taking measurements after each conditioning procedure. All parameters were kept constant for conditioning as in Section 2.3.2, although the position of the conditioning ring was changed to affect the plate shape.





Figure 7,8: Plot of the plate surface following conditioning for 35 minutes. Figure 7 (above) shows the plate surface is beginning to even out, although the edges and center are still higher than the middle portion of the plate area. Adjustment of the conditioning ring is needed to correct the problem. Figure 8 (below) shows the plate after conditioning in the new position for 37 minutes. The plate is now within specification.



4.0: Results

It has been demonstrated that proper plate conditioning is important to maintaining good specimens when flatness is of utmost importance. Cast iron lapping plates, such as the LP 920M, require conditioning almost after every process run to make sure that the plates are maintained. Conditioning during lapping will undoubtedly minimize the amount of time needed for conditioning the plate individually, and will enhance the resulting specimens.

1. Plate conditioning is required when using cast iron lapping plates with free abrasives.
2. The existing conditioning ring may not be the most efficient method for conditioning as it seems to have a negligible effect on the plate surface after a considerable amount of lapping time.
3. Conditioning the plates using a heavier, thicker conditioning ring reduced conditioning time and improved the plate surface flatness.
4. Following the use of lapping and polishing fixtures it appears to be necessary to re-condition the lapping plate to maintain good flatness.
5. A repeatable method for plate flatness inspection was developed.

